

# Denture appointment armamentarium

	Traditional	Additional if implant supported/retained
<b>Diagnostic impressions</b>	Alginate Fox plane (sterile) Ruler	Diagnostic T (sterile)
<b>Final impression</b>	Custom tray(s) Tray adhesive Glass well Brushes Heavy & medium body PVS material #21 laboratory blade & handle Thompson markers (if maxillary) Iowa wax (if maxillary)	Implant kit Impression posts Light body PVS material Floss Pattern resin Glass wells x2 Brushes Radiograph
<b>Jaw relations</b>	Wax rims on trial bases from lab Flame & butane refill Fox plane (sterile) Facebow transfer (earpiece, jig, bitefork (sterile), facebow ruler) Fast set PVS Buffalo knife Ruler Tooth book & shade guide	Implant kit  <b>If verification jig appt also:</b> Radiograph Pattern resin Laboratory wells x2 Brushes
<b>Wax try-in</b>	Set-up and mounted models from lab Horseshoe articulating paper Fox plane (sterile) Ruler Hand mirror	Implant kit  <b>If verification jig appt also:</b> Radiograph Pattern resin Laboratory wells x2 Brushes
<b>Delivery</b>	Dentures Pressure indicating paste & brush Horseshoe articulating paper Straight nose cone and acrylic bur (sterile) Hand mirror	Implant kit Radiograph Digital photos
<b>Post-op denture adjustment</b>	Pressure indicating paste & brush Horseshoe articulating paper Straight nose cone and acrylic bur	